

			FrameStar® 2-Component PCR Plates							FrameStar® PCR Plates, RA Plate, FrameStrips®								Standard PCR Plates and Strips																			
Skirt			skirted		semi-skirted					non skirted		Break-A-Way		Break-2-Ways		skirted	strips	skirted		semi-skirted			non	strips			strips										
Number of wells			384	384	96	192	96	96	96	96	96	96	96	96	96	96	8	384	384	96	96	96	96	96	8	8	8	8	8	4							
Standard/Low Profile			/	/	L	/	L	L	S	S	S	S	L	L	L	S	/	/	L	S	L	S	S	S	L	S	L	S	/								
Details see PCR Consumables Overview Brochure page			8	8	9	8	9	10	10	10	11	11	11	13	13	13	13	15	16	18	18	18	19	19	19	20/22	20	17	21	21	21						
Product code ¹ /Short description			4ti-0384 FrameStar® 384	4ti-0380 FrameStar® 384 Roche	4ti-0960 FrameStar® 96	4ti-0192 FrameStar® 192	4ti-0950 & 4ti-0954 FrameStar® 96 Roche	4ti-0910 FrameStar® 96 FastPlate	4ti-0730 FrameStar® 96 ABI®	4ti-0770 FrameStar® 96 ABI®	4ti-0900 FrameStar® 96	4ti-0710 FrameStar® 96	4ti-0720 FrameStar® 96	4ti-1000 FrameStar® Break-A-Way	4ti-1200 FrameStar® Break-A-Way	4ti-1300 FrameStar® Break-2-Ways	4ti-1400 FrameStar® Break-2-Ways	4ti-0960/RA Random Access 96	4ti-0785 FrameStrip® 8 Well	4ti-1384 384 Well	4ti-1381 384 Well Roche	4ti-0740 96 Well	4ti-0760 96 Well	4ti-0955 96 Well Roche	4ti-0735 96 Well ABI®	4ti-0750 - 96 Well 4ti-0750/TA- dividable	4ti-0781 8 Well	4ti-0753 Vari-Strips™ 8 Well ²	4ti-0792 8 Well, Attached Flat Caps	4ti-0793 8 Well, Attached Flat Caps	4ti-0794 8 Well, Attached Domed C.	4ti-0796 4 Well Rotor-Gene®					
ABI® / LIFE TECHNOLOGIES / THERMO FISHER SCIENTIFIC																																					
Thermal Cyclers	96 well standard block	Veriti, Proflex, Simpliamp							✓	Ⓢ				✓																							
	96 well FAST block	GeneAmp® 2700 / 2720 / 9600 / 9700							✓	Ⓢ				✓																							
	384 well block	GeneAmp® 9700, Veriti, Proflex, Multiblock system	Ⓢ			•															✓																
qPCR Cyclers	96 well standard block	7000, 7300, 7500, 7700, 7900 HT							✓	Ⓢ				✓																							
		QuantStudio™ 3 / 5 / 6 / 7 / 12K, ViiA7™							✓	Ⓢ				✓																							
		StepOne												✓																							
	96 well FAST block	7500 FAST, 7900 HT FAST							Ⓢ																												
Sequencers	384 well block	QuantStudio™ 5 / 6 / 7 / 12K, ViiA7™, 7900 HT FAST	Ⓢ						Ⓢ																												
	96 well block	3100, 3130XL, 3500, 3500XL, 3730, 3730XL							Ⓢ	✓																											
	384 well block	3100, 3130XL, 3500, 3500XL, 3730, 3730XL	Ⓢ																																		
AGILENT / STRATAGENE																																					
Thermal Cyclers	96 well block	Surecycler 8800																																			
	384 well block	SureCycler 8800	Ⓢ			•																															
	96 well block	Robocycler Gradient			Ⓢ																																
qPCR Cyclers	96 well block	AriaMx			Ⓢ																																
		Mx3000P™, Mx3005P™								Ⓢ	✓																										
		Mx4000™							✓	Ⓢ																											
ANALYTIK JENA / BIOMETRA																																					
Thermal Cyclers	Strips only	TRIO, Tpersonal, T3 Thermocycler																																			
	96 well block	Flexcycler2, T1 Thermocycler, Tgradient, Tone, Tadvanced, TProfessional (Standard/Basic) Gradient/XL			Ⓢ																																
	384 well block	Trobot 96, SpeedCycler2 (SP, SPR)			Ⓢ																																
qPCR Cyclers	96 well block	Flexcycler2, T1 Thermocycler, Tadvanced, TProfessional, Trobot 96	Ⓢ			•																															
	384 well block	qTOWER ³ / G / touch, Toptical																																			
	384 well block	qTOWER ³ 84 / 84G	Ⓢ			•																															
BIOER TECHNOLOGIES																																					
Thermal Cyclers	Strips only	GeneQ																																			
	96 well block	Gene Touch 96																																			
	384 well block	Gene Touch 384	Ⓢ			•																															
BIO-RAD																																					
Thermal Cyclers	Strips only	Genecycler																																			
	96 well block	C1000 Touch, S1000			Ⓢ																																
	384 well block	iCycler™, MyCycler™, T100																																			
qPCR Cyclers	96 well block	C1000 Touch, S1000	Ⓢ			•																															
	384 well block	CFX96 Touch, CFX96 Touch Deep Well, CFX connect			Ⓢ																																
	384 well block	MyiQ™, iCycler™ IQ / IQ 4 / IQ 5																																			
	384 well block	CFX384 Touch	Ⓢ			•																															
BIO-RAD MJ RESEARCH																																					
Thermal Cyclers	Strips only	Mini Gradient																																			
	96 well block	Personal																																			
	384 well block	PTC100™ / 200™ / 220™ / 221™ / 225™ / 240™			Ⓢ																																
	384 well block	PTC200™ / 220™ / 221™ / 225™ / 240™	Ⓢ			•																															

¹ Short product code shown only (without details on frame and well colour), please refer to the corresponding product page for ordering details of all variations available.
² For compatibility information of the Vari-Plates™ please refer to the table entry for the respective Vari-Plate™ frame.
 Ⓢ Recommended option ✓ Compatible • Should be compatible, please check with your specific instrument/block TA Compatible with separated Tear-A-Way™ Plates, only

The individual instrument-plate compatibility information is for guidance only. Samples of all plates are available before purchase to ensure compatibility. Please check installed heat block, refer to instruments manual for details or contact our technical support.
 4titude® recognises that designated trademarks and brands of the Instrument Compatibility Table are the property of their respective owners.

		FrameStar® 2-Component PCR Plates								FrameStar® PCR Plates, RA Plate, FrameStrips®								Standard PCR Plates and Strips															
Skirt		skirted				semi-skirted				non skirted		Break-A-Way		Break-2-Ways		skirted	strips	skirted		semi-skirted		non	strips				strips						
Number of wells		384	384	96	192	96	96	96	96	96	96	96	96	96	96	96	8	384	384	96	96	96	96	8	8	8	8	8	4				
Standard/Low Profile		/	/	L	/	L	L	S	S	S	S	L	S	L	L	S	/	/	L	S	L	S	S	S	L	S	L	S	/				
Details see PCR Consumables Overview Brochure page		8	8	9	8	9	10	10	10	11	11	11	13	13	13	13	15	16	18	18	18	19	19	19	20/22	20	17	21	21	21			
Product code ¹ /Short description		4ti-0384 FrameStar® 384	4ti-0380 FrameStar® 384 Roche	4ti-0960 FrameStar® 96	4ti-0192 FrameStar® 192	4ti-0950 & 4ti-0954 FrameStar® 96 Roche	4ti-0910 FrameStar® 96 FastPlate	4ti-0730 FrameStar® 96 ABI®	4ti-0770 FrameStar® 96 ABI®	4ti-0900 FrameStar® 96	4ti-0710 FrameStar® 96	4ti-0720 FrameStar® 96	4ti-1000 FrameStar® Break-A-Way	4ti-1200 FrameStar® Break-A-Way	4ti-1300 FrameStar® Break-2-Ways	4ti-1400 FrameStar® Break-2-Ways	4ti-0960/RA Random Access 96	4ti-0785 FrameStrip® 8 Well	4ti-1384 384 Well	4ti-1381 384 Well Roche	4ti-0740 96 Well	4ti-0760 96 Well	4ti-0955 96 Well Roche	4ti-0735 96 Well ABI®	4ti-0750 - 96 Well 4ti-0750/TA- dividable	4ti-0781 8 Well	4ti-0753 Vari-Strips™ 8 Well ²	4ti-0792 8 Well, Attached Flat Caps	4ti-0793 8 Well, Attached Flat Caps	4ti-0794 8 Well, Attached Domed C.	4ti-0796 4 Well Rotor-Gene®		
BIO-RAD MJ RESEARCH <i>continued</i>																																	
qPCR Cylers	Strips only			4							✓	✓	✓	✓	✓	✓	✓																
	96 well block			4					✓	✓															✓	✓	4	✓	✓				
CORBETT RESEARCH																																	
Thermal Cylers	96 well block								4	✓		✓																					
	384 well block	4			•														✓														
qPCR Cycl.	Strips only																															4	
EPPENDORF																																	
Thermal Cylers	96 well block			4							✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	✓		
	384 well block	4			•							4	✓	✓	✓	✓	✓	✓							✓	✓	✓	✓	✓	✓	✓		
qPCR Cycl.	96 well block			4							✓	✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	✓		
GE HEALTHCARE / AMERSHAM																																	
Sequencers	96 well block			4								✓		✓		✓					✓												
	384 well block	4			•														✓														
PEQLAB / VWR																																	
Thermal Cylers	Strips only																																
	96 well block			4							✓	✓	✓	✓	✓	✓	✓	✓							TA	✓	✓	✓	✓	✓	✓		
	384 well block	4			✓														✓		✓	✓			✓	✓	✓	✓	✓	✓	✓		
ROCHE																																	
qPCR Cylers	96 well block					4																✓					✓						
	384 well block		4																	✓													
	Strips only																									4							
SENSOQUEST																																	
Thermal Cylers	96 well block			✓			✓	✓	✓	✓		✓	✓	✓	✓	✓	✓	✓			✓	✓			✓	✓	✓	✓	✓	✓	✓		
	192 well block					✓																											
	384 well block	4			•															✓													
TAKARA																																	
Thermal Cycl.	96 well block								4	✓		✓		✓			✓				✓				✓	✓		✓	✓	✓	✓		
TECHNE																																	
Thermal Cylers	Strips only											4		✓			✓								TA	✓		✓		✓			
	96 well block										4	✓	✓	✓	✓	✓	✓	✓				✓			✓	✓	✓	✓	✓	✓	✓		
	384 well block	4			•														✓														
qPCR Cycl.	96 well block	4		4	✓														✓		✓						✓		✓				

¹ Short product code shown only (without details on frame and well colour), please refer to the corresponding product page for ordering details of all variations available.
² For compatibility information of the Vari-Plates™ please refer to the table entry for the respective Vari-Plate™ frame.

4 Recommended option ✓ Compatible • Should be compatible, please check with your specific instrument/block TA Compatible with separated Tear-A-Way™ Plates, only

The individual instrument-plate compatibility information is for guidance only. Samples of all plates are available before purchase to ensure compatibility. Please check installed heat block, refer to instruments manual for details or contact our technical support.
 4titude® recognises that designated trademarks and brands of the Instrument Compatibility Table are the property of their respective owners.